

VIBRATING TABLE

Code : B125



- Designed from robust steel and has a very steady frame
- Used to compact the concrete specimens in the laboratory
- Equipped with a vibrating motor
- Supplied with a clamp assembly for moulds.

STANDARDS

ASTM C192 • EN 12390-2

ORDERING INFORMATION

Item Name	Item Code
CLAMPING APPARATUS for VIBRATING TABLE (Small)	B125P001H

CLAMPING APPARATUS for VIBRATING TABLE (Large)	B125P002H
VIBRATING TABLE (Small)	B125X001U
VIBRATING TABLE (Large)	B125X002U
VIBRATING TABLE (Large) - 380V	B125X002V