

CONCRETE > SETTING TIME BY PENETRATION

CONCRETE DIAL PENETROMETER

Code : B308



Used to determine the initial set of concrete mortar

STANDARDS

ASTM C403 • AASHTO T197 • UNI 7123

TECHNICAL SPECIFICATIONS

- Capacity: 50 kgf (700 lbs/sq.in = 700 psi)
- Dial with scale readable in both (kgf/cm²) and (lbf/sq.in)
- Plunger's surface area: 32.30 mm² (~ 1/20 sq.in)

ORDERING INFORMATION

Item Name	Item Code
CONCRETE INITIAL SET PENETROMETER - Dial Model	B308X00XH