

CEMENT COMPRESSION and FLEXURAL TESTING MACHINE (with Loadcell)

Code : C005



will be added soon ...

TECHNICAL SPECIFICATIONS

will be added soon ...

ORDERING INFORMATION

Item Name	Item Code
HIGH-CAPACITY CEMENT COMPRESSION and FLEXURE TESTER (15/250 kN, with Loadcell)	C005A25XU
HIGH-CAPACITY CEMENT COMPRESSION and FLEXURE TESTER (25/600 kN, with Loadcell)	C005A60XU