

CEMENT > WORKABILITY & CONSISTENCY

CEMENT FLOW TABLE (Motorized)

Code : C021M



will be added soon ...

TECHNICAL SPECIFICATIONS

will be added soon ...

ORDERING INFORMATION

Item Name	Item Code
CEMENT FLOW TABLE - ASTM	C021H00AH
CEMENT FLOW TABLE - EN	C021H00EH
MOTORIZED CEMENT FLOW TABLE - ASTM [60 Hz]	C021M00AK
MOTORIZED CEMENT FLOW TABLE - ASTM [50 Hz]	C021M00AT
MOTORIZED CEMENT FLOW TABLE - EN [60 Hz]	C021M00EK

MOTORIZED CEMENT FLOW TABLE - EN [50 Hz]	C021M00ET
BRASS CONICAL MOULD - ASTM	C021P001H
BRASS CONICAL MOULD - EN	C021P002H
TAMPER - ASTM	C021P003H
TAMPER - EN	C021P006H