

CEMENT MIXER (Automatic, Programmable)

Code : C050A



will be added soon ...

TECHNICAL SPECIFICATIONS

will be added soon ...

ORDERING INFORMATION

Item Name	Item Code
AUTOMATIC PROGRAMMABLE CEMENT MIXER - 5 lt	C050A00XU
MORTAR MIXER - 5 lt	C050M00XU
MORTAR MIXER - 10 lt	C050M10XU
BEATER for MORTAR MIXER - 5 lt	C050P001H
STEEL BOWL for MORTAR MIXER - 5 lt	C050P002H

BEATER (Stainless Steel) for MORTAR MIXER - 5 lt

C050P003H

ELECTRONIC BOARD (PCB) for C050A

C050P004H