

ELECTROMAGNETIC SIEVE SHAKER (Analog)

Code : G041M



will be added soon ...

TECHNICAL SPECIFICATIONS

will be added soon ...

ORDERING INFORMATION

Item Name	Item Code
ELECTROMAGNETIC SIEVE SHAKER - Ø 200 mm	G041A20XU
ELECTROMAGNETIC SIEVE SHAKER - Ø 300 mm	G041A30XU
ELECTROMAGNETIC SIEVE SHAKER - Ø 450 mm	G041A45XU
ELECTROMAGNETIC SIEVE SHAKER, ANALOG - Ø 200 mm	G041M20XU
ELECTROMAGNETIC SIEVE SHAKER, ANALOG - Ø 300 mm	G041M30XU

ELECTROMAGNETIC SIEVE SHAKER, ANALOG - Ø 450 mm	G041M45XU
NOISE-REDUCTION and SAFETY CABINET for SIEVE SHAKER Ø 200 mm	G041P001H
NOISE-REDUCTION and SAFETY CABINET for SIEVE SHAKER Ø 300 mm	G041P002H
NOISE-REDUCTION and SAFETY CABINET for SIEVE SHAKER Ø 450 mm	G041P003H
SPARE FIBER PLATES SET for ELECTROMAGNETIC SIEVE SHAKERS	G041P004H
SPARE KNOB for SIEVE SHAKER	G041P005H
SPARE DISPLAY PANEL (4.3") for G041A	G041P006H
ELECTRONIC BOARD (PCB) for G041	G041P007H
CLAMP PLATE for G041A45XU	G041P008H