

VACUUM PUMP

Code : G702



Used to vacuum specimens, remove air voids and speed up deairing process.

ORDERING INFORMATION

Item Name	Item Code
HOSE for VACUUM PUMP	G702P001H
GAUGE for VACUUM PUMP	G702P002H
REGULATOR for VACUUM PUMP	G702P003H
VACUUM PUMP - 51 lt/min	G702X051U
VACUUM PUMP (Oil-Free) - 99 mBar / 62 lt/min	G702X062U
VACUUM PUMP - 100 lt/min	G702X100U
VACUUM PUMP - 128 lt/min (Dual Stage)	G702X128U

VACUUM PUMP - 170 lt/min (Dual Stage)

G702X170U