

GLASS THERMOMETER

Code : G780



- Thermometers in engineering testing laboratories are used to measure the temperature of samples, materials, or equipment during experiments. They help ensure accurate temperature control and monitoring, which is critical for processes that are temperature-sensitive, such as curing, heating, or material testing.

ORDERING INFORMATION

Item Name	Item Code
GLASS THERMOMETER - 100°C	G780X100H
GLASS THERMOMETER - 200°C	G780X200H
GLASS THERMOMETER - 250°C	G780X250H
GLASS THERMOMETER - 360°C	G780X360H