

MULTI-DECK GRADING MACHINE

Code : M042



- Used specially in mining labs where large quantities of sieving is required.
- Two sieves with decks, each of 300 x 500 mm size, and the openings are 10 mm and 5 mm unless otherwise specified during the order.
- Plastic boxes, as collectors, are supplied for each deck.

TECHNICAL SPECIFICATIONS

- Equipped with pan at the bottom.
- Adjustable sieving angle.
- 380 V / 50 Hz

ORDERING INFORMATION

Item Name	Item Code
SIEVE MESH, 400 x 800 mm, 0.250 mm (No.60)	M0420250Y
SIEVE MESH, 400 x 800 mm, 0.500 mm (No.35)	M0420500Y
SIEVE MESH, 400 x 800 mm, 1.00 mm (No.18)	M0421000Y
MULTI-DECK GRADING MACHINE	M042X00XU