

# ABRASION and CRUMBLING INDEX TESTING MACHINE

Code : M191



- Used to determine the resistance to abrasion and crumbling of the iron ores, such as pellet and sinter form
- The drum rotates 200 revolutions with 15 kg of sample
- Supplied with 2 test sieves
- Collector at the bottom to collect shredded materials

## **TECHNICAL SPECIFICATIONS**

- The drum rotates 200 revolutions with 15 kg of sample
- Collector at the bottom to collect shredded materials
- 220 - 240 V / 50 - 60 Hz

**ORDERING INFORMATION**

| Item Name                           | Item Code |
|-------------------------------------|-----------|
| ABRASION and CRUMBLING INDEX TESTER | M191X00XU |