

HARDGROVE GRINDABILITY INDEX TESTING MACHINE

Code : M195



- Used to determine the HGI (Hardgrove Grindability Index) of coal
- Supplied with abrasive spheres (8 ea)
- Weight set is supplied with the machine

TECHNICAL SPECIFICATIONS

- Digital numerator for setting the revolution of the machine and automatically stopping after desired revolution
- 220 - 240 V / 50 - 60 Hz

ORDERING INFORMATION

Item Name	Item Code
HARDGROVE GRINDABILITY INDEX TESTER	M195X00XU