


SOIL > ATTERBERG LIMITS

LIQUID LIMIT DEVICE (Semi-Automatic, Cone Penetrometer Method)

Code : T090S



will be added soon ...

TECHNICAL SPECIFICATIONS

will be added soon ...

ORDERING INFORMATION

Item Name	Item Code
DIGITAL CONE PENETROMETER	T090H00XB
SPARE DIGITAL INDICATOR for PENETROMETER	T090P001H
PENETRATION CONE - 30° x 80 g - BS	T090P002H

PENETRATION CONE - 60° x 60 g - EN	T090P003H
PENETRATION CONE - 30° x 80 g - EN	T090P004H
PENETRATION CONE - 30° x 100 g - EN	T090P005H
PENETRATION CONE - 30° x 400 g - EN	T090P006H
CONE TEST GAUGE - 30°	T090P007H
CONE TEST GAUGE - 60°	T090P008H
SLIDING ROD for T090	T090P009H
SEMI-AUTOMATIC DIGITAL CONE PENETROMETER	T090S00XU