

SOIL > SHEAR STRENGTH

DIRECT SHEAR TESTING MACHINE (Large, Conventional)

Code : T220 (Large)



will be added soon ...

TECHNICAL SPECIFICATIONS

will be added soon ...

ORDERING INFORMATION

Item Name	Item Code
DIRECT SHEAR TESTER (Conventional Model)	T220A10XU
LARGE DIRECT SHEAR TESTER - 300 mm	T220A30XU
LARGE DIRECT SHEAR TESTER - 500 mm	T220A50XU

SHEAR BOX ASSEMBLY - Ø 60 mm	T220P001H
SHEAR BOX ASSEMBLY - Ø 100 mm	T220P002H
SHEAR BOX ASSEMBLY - 60 x 60 mm	T220P003H
SHEAR BOX ASSEMBLY - 100 x 100 mm	T220P004H
LOADING TOP CAP - Ø 60 mm(for DIRECT SHEAR)	T220P005H
POROUS DISCS - Ø 60 mm (for DIRECT SHEAR)	T220P006H
POROUS STONES - Ø 60 mm (for DIRECT SHEAR)	T220P007H
CUTTING EDGE - Ø 60 mm (for DIRECT SHEAR)	T220P008H
EXTRUSION DOLLY - Ø 60 mm (for DIRECT SHEAR)	T220P009H
LOADING TOP CAP - 60x60 mm(for DIRECT SHEAR)	T220P010H
POROUS DISCS - 60x60 mm (for DIRECT SHEAR)	T220P011H
POROUS STONES - 60x60 mm (for DIRECT SHEAR)	T220P012H
CUTTING EDGE - 60x60 mm (for DIRECT SHEAR)	T220P013H
EXTRUSION DOLLY - 60x60 mm (for DIRECT SHEAR)	T220P014H