

STANDALONE DE-AIRING SYSTEM

Code : T235



- Used to remove dissolved air from water before it is used in soil mechanics testing, such as triaxial, permeability, and related geotechnical test methods

AVAILABLE MODELS

Standalone De-Airing System:

- Capacity: 20 lt
- Equipped with:
 - Magnetic stirrer
 - Front valve panel
 - Tube inlet and outlet
 - Vacuum Pump

- Tubes outer diameter: Ø 6 mm
- 220 - 240 V / 50 - 60 Hz

Standalone De-Airing Tank:

- Capacity: 7 lt
- Wall-Mounted
- Equipped with 3 valves:
 - Water Inlet/Outlet Valve (at the bottom)
 - Vacuum Pump Inlet Valve (at the top)
 - Pressure Release Valve (at the top)
- Tubes outer diameter: Ø 6 mm
- Tubes and vacuum pump should be ordered separately

ORDERING INFORMATION

Item Name	Item Code
STANDALONE DE-AIRING TANK - 7 lt	T235H007H
STANDALONE DE-AIRING SYSTEM - 20 lt	T235X020U

OTHER PHOTOS

