

SAMPLE EXTRUDER (Universal, Horizontal)

Code : U020 (H)



- This tool is designed for extracting samples from molds, tubes, Shelby tubes, and cutting rings.
- It features a hydraulic system for quick and effortless extraction with the press of a single button.
- With an impressive 900 mm ram travel, it ensures efficient operation, complemented by speed control for precise handling.

STANDARDS

AASHTO T134 • AASHTO T180 • AASHTO T193 • AASHTO T245 • ASTM D1557 • ASTM D1559 • ASTM D1587-15 • ASTM D1883 • ASTM D5102 • ASTM D558 • ASTM D560 • ASTM D698 • BS 1377:4 • BS 598:107 • EN 12697-30 • EN 13286-2 • EN 13286-47

EQUIPPED WITH

- Hydraulic System

SUPPLIED WITH

- Ø 100 mm adaptors set
- Ø 150 mm adaptors set

ORDERING INFORMATION

Item Name	Item Code
UNIVERSAL SAMPLE EXTRUDER - 100 mm	U020H100H
UNIVERSAL SAMPLE EXTRUDER - 300 mm	U020H300H
UNIVERSAL SAMPLE EXTRUDER - 450 mm	U020H450H
UNIVERSAL SAMPLE EXTRUDER - 600 mm	U020H600H
MOTORIZED SAMPLE EXTRUDER (Desktop Version) - 100 mm	U020M100U
MOTORIZED SAMPLE EXTRUDER - 300 mm	U020M300U
MOTORIZED SAMPLE EXTRUDER - 650 mm	U020M650U
MOTORIZED SAMPLE EXTRUDER (Horizontal) - 900 mm	U020M900U
ADAPTOR for [U 020] - Ø 38 mm	U020P001H
ADAPTOR for [U 020] - Ø 50 mm	U020P002H
ADAPTOR for [U 020] - Ø 70 mm	U020P003H
ADAPTOR for [U 020] - Ø 106 mm	U020P004H
JACK (30 kN) for U020	U020P005H
ADAPTORS SET for HORIZONTAL EXTRUDER (Ø 50.80 / 48.50 / 43.50 mm)	U020PH01H
ADAPTORS SET for HORIZONTAL EXTRUDER (Ø 63.50 / 61.00 / 56.00 mm)	U020PH02H
ADAPTORS SET for HORIZONTAL EXTRUDER (Ø 75.00 / 72.00 / 67.00 mm)	U020PH03H

ADAPTORS SET for HORIZONTAL EXTRUDER (Ø 83.00 / 79.50 / 74.50 mm)	U020PH04H
ADAPTORS SET for HORIZONTAL EXTRUDER (Ø 88.90 / 85.00 / 80.00 mm)	U020PH05H
ADAPTORS SET for HORIZONTAL EXTRUDER (Ø 100.00 / 96.00 / 91.00 mm)	U020PH06H
ADAPTORS SET for HORIZONTAL EXTRUDER (Ø 101.60 / 97.50 / 92.50 mm)	U020PH07H
ADAPTORS SET for HORIZONTAL EXTRUDER (Ø 114.30 / 109.00 / 104.00 mm)	U020PH08H
ADAPTORS SET for HORIZONTAL EXTRUDER (Ø 125.00 / 119.00 / 114.00 mm)	U020PH09H
ADAPTORS SET for HORIZONTAL EXTRUDER (Ø 106.00 / 100.00 / 95.00 mm)	U020PH10H
ADAPTORS SET for HORIZONTAL EXTRUDER (Ø 156.00 / 150.00 / 145.00 mm)	U020PH11H