

TEMPERATURE MONITORING and RECORDING SYSTEM

Code : U030



will be added soon ...

TECHNICAL SPECIFICATIONS

will be added soon ...

ORDERING INFORMATION

Item Name	Item Code
TEMPERATURE SENSOR (PT100)	U030P001H
TEMPERATURE SENSOR (Water-Proof, DS18B20)	U030P002H
TEMPERATURE & HUMIDITY SENSOR	U030P003H

TEMPERATURE SENSOR (K-Type Thermocouple)	U030P004H
TEMPERATURE SENSOR (LAE NTC SN4B20P1)	U030P005H
TEMPERATURE MONITORING and RECORDING SYSTEM	U030X00XU