

LE CHATELIER FLASK

Code : W019



- Le Chatelier Flask is used to determine the relative density (the specific gravity) of hydraulic cement and lime.
- The measurement is done by taking an empty flask, filling it with cement then adding polar liquid and measure the weight of the flask at each step and using it in the related formula.

TECHNICAL SPECIFICATIONS

- Graduation (0-1) & (18-24) ml. with 0.05 ml. accuracy.

ORDERING INFORMATION

Item Name	Item Code
LE CHATELIER FLASK - 250 ml	W019G250H