

PORCELAIN MORTAR and RUBBER-HEADED PESTLE

Code : W027



- A porcelain mortar and rubber-headed pestle in engineering testing laboratories are used for grinding and mixing materials without causing contamination or damage. The rubber-headed pestle is gentle on the mortar's surface and is ideal for handling delicate or fragile samples, making it useful for preparing samples in chemical, biological, or material testing.

ORDERING INFORMATION

Item Name	Item Code
PORCELAIN MORTAR and RUBBER-HEADED PESTLE	W027P001H