

PLASTIC BEAKER

Code : W034



- A plastic beaker in engineering testing laboratories is used for holding, mixing, and heating liquids or chemicals. It is lightweight, durable, and resistant to breaking, making it ideal for everyday use in various laboratory procedures where glass beakers may be too fragile or impractical.

ORDERING INFORMATION

Item Name	Item Code
PLASTIC BEAKER - 1000 ml	W034P01KH
PLASTIC BEAKER - 25 ml	W034P025H
PLASTIC BEAKER - 2000 ml	W034P02KH
PLASTIC BEAKER - 50 ml	W034P050H
PLASTIC BEAKER - 100 ml	W034P100H

PLASTIC BEAKER - 250 ml

W034P250H

PLASTIC BEAKER - 500 ml

W034P500H