

PLASTIC FUNNEL

Code : W035



- A plastic funnel in engineering testing laboratories is used for transferring liquids or powders into containers with narrow openings. It is lightweight, durable, and resistant to breaking, making it ideal for handling various substances in laboratory processes while preventing spills and ensuring controlled flow.

ORDERING INFORMATION

Item Name	Item Code
PLASTIC FUNNEL - 50 mm	W035P050H
PLASTIC FUNNEL - 70 mm	W035P070H
PLASTIC FUNNEL - 90 mm	W035P090H
PLASTIC FUNNEL - 100 mm	W035P100H
PLASTIC FUNNEL - 120 mm	W035P120H

PLASTIC FUNNEL - 150 mm

W035P150H