

GLASS PLATE

Code : W153



- A glass plate in engineering testing laboratories is used for various purposes, including surface preparation, holding samples, or as a work surface for experiments. It provides a smooth, non-reactive surface, making it ideal for activities like chemical reactions, chromatography, or for observing materials under a microscope.

ORDERING INFORMATION

Item Name	Item Code
GLASS PLATE - 50 x 50 mm	W153G050H
GLASS PLATE - 300 x 300 mm	W153G300H
GLASS PLATE - 500 x 500 mm	W153G500H