

ASPHALT > DURIEZ

DURIEZ TEST SET

Code : A003



will be added soon ...

STANDARDS

will be added soon ...

TECHNICAL SPECIFICATIONS

will be added soon ...

ORDERING INFORMATION

Item Name	Item Code
-----------	-----------

DURIEZ TESTING MOULD - Ø 80 mm	A003P001H
PENETRATION PISTON - Ø 80 mm	A003P002H
UPPER and LOWER PISTON - Ø 80 mm	A003P003H
TEMP. SUPPORT - Ø 80 mm	A003P004H
DEMOULDING CYLINDRICAL CONTAINER - Ø 80 mm	A003P005H
DURIEZ TESTING MOULD - Ø 120 mm	A003P006H
PENETRATION PISTON - Ø 120 mm	A003P007H
UPPER and LOWER PISTON - Ø 120 mm	A003P008H
TEMP. SUPPORT - Ø 120 mm	A003P009H
DEMOULDING CYLINDRICAL CONTAINER - Ø 120 mm	A003P010H
DURIEZ TEST SET - Ø 80 mm	A003X08FH
DURIEZ TEST SET - Ø 120 mm	A003X12FH

OTHER PHOTOS

