

ASPHALT > MARSHALL STABILITY

MARSHALL MOULD

Code : A010



- Used in preparing and compacting asphalt specimens. The prepared specimen is then used in Marshall Stability tests to determine the strength and flow of that specific Marshall Mix.
- The mould is made of a heavy-duty steel and protected against corrosion.
- The mould is supplied with base plate and collar.

STANDARDS

ASTM D6926 • EN 12697-10 • AASHTO t245

TECHNICAL SPECIFICATIONS

- Made of Galvanized Steel
- Diameter:

- Ø 4" B-010/4
- Ø 6" B-010/6

ORDERING INFORMATION

Item Name	Item Code
COLLAR for MARSHALL MOULD - 4"	A010P001H
MOULD BODY for MARSHALL MOULD - 4"	A010P002H
BASE PLATE for MARSHALL MOULD - 4"	A010P003H
COLLAR for MARSHALL MOULD - 6"	A010P004H
MOULD BODY for MARSHALL MOULD - 6"	A010P005H
BASE PLATE for MARSHALL MOULD - 6"	A010P006H
FILTER PAPER for MARSHALL MOULD - 4"	A010P007H
FILTER PAPER for MARSHALL MOULD - 6"	A010P008H
MARSHALL MOULD - 4"	A010X04AH
MARSHALL MOULD - 6"	A010X06AH