

ASPHALT > MARSHALL STABILITY

MARSHALL COMPACTOR (Manual)

Code : A016



- Used to determine the resistance of an asphalt specimen to plastic flow.
- The test is performed using a manually dropped hammer.

STANDARDS

ASTM D6926

TECHNICAL SPECIFICATIONS

Composed of:

- Compaction Hammer
- Wooden Pedestal
- Steel Plate

- Mould Holder
- Hammer Guide

ORDERING INFORMATION

Item Name	Item Code
MARSHALL COMPACTOR, HAND-OPERATED - ASTM	A016H04AH

OTHER PHOTOS

