

## REFLUX EXTRACTOR

Code : A035



- Used to determine the percentage of bitumen in the asphalt pavement mixtures.
- Equipped with a glass cylinder, two wire mesh cones and brass condenser.
- Filter Papers should be ordered separately

### **STANDARDS**

ASTM D2172 • AASHTO T-164

### **TECHNICAL SPECIFICATIONS**

| Code               | A-035/1   | A-035/4   |
|--------------------|-----------|-----------|
| Capacity           | 1000 g    | 4000 g    |
| Glass Jar Cylinder | A-035/1/G | A-035/4/G |

Wire Mesh Cone (Set of 2)

Thermometer

Hotplate

GL-11/S1

GL-11/S1

Brass Condenser

### ORDERING INFORMATION

| Item Name                                            | Item Code |
|------------------------------------------------------|-----------|
| CONDENSER for A-035/1                                | A035P001H |
| REFLUX BITUMEN EXTRACTOR GLASS - 1000 g              | A035P002H |
| MESH CONES SET for REFLUX BITUMEN EXTRACTOR - 1000 g | A035P003H |
| CONDENSER for A-035/4                                | A035P004H |
| REFLUX BITUMEN EXTRACTOR GLASS - 4000 g              | A035P005H |
| MESH CONES SET for REFLUX BITUMEN EXTRACTOR - 4000 g | A035P006H |
| FILTER PAPER - Ø 40 cm                               | A035P007H |
| REFLUX BITUMEN EXTRACTOR - 1000 g                    | A035X10XU |
| REFLUX BITUMEN EXTRACTOR - 4000 g                    | A035X40XU |