

ASPHALT > CORING

CORE DRILLING MACHINE (Trailer Mounted)

Code : A065



- Used to take core samples from asphalt.
- Installed on a trailer for easy transportation.
- Equipped with a petrol engine, four stabilizing feet, 100 lt water tank.
- Core Bits for the machine should be ordered separately

STANDARDS

ASTM D 5361 • EN 12697-27

TECHNICAL SPECIFICATIONS

- Trailer Mounted

- Petrol Engine
- Stabilizing Feet (4 ea)
- Water Tank (100 Lt)

ORDERING INFORMATION

Item Name	Item Code
CORE BIT for ASPHALT PAVEMENT - Ø 50 mm x 400 mm	A065P001H
PETROL ENGINE for CORE DRILLING MACHINE	A065P001P
CORE BIT for ASPHALT PAVEMENT - Ø 75 mm x 400 mm	A065P002H
CORE BIT for ASPHALT PAVEMENT - Ø 100 mm x 400 mm	A065P003H
CORE BIT for ASPHALT PAVEMENT - Ø 150 mm x 400 mm	A065P004H
CORE BIT for ASPHALT PAVEMENT - Ø 200 mm x 400 mm	A065P005H
FRAME for CORE DRILLING MACHINE (A065)	A065P006H
WATER TANK for CORE DRILLING MACHINE (A065)	A065P007H
CORE DRILLING MACHINE - 4 HP (Trailer Mounted)	A065X00XP

OTHER PHOTOS

