

BITUMEN PENETROMETER (Automatic)

Code : A090A



will be added soon ...

TECHNICAL SPECIFICATIONS

will be added soon ...

ORDERING INFORMATION

Item Name	Item Code
BITUMEN PENETROMETER - Automatic	A090A00XU
BITUMEN PENETROMETER - Digital	A090H00XB
NEEDLE for BITUMEN PENETROMETER	A090P001H
TRANSFER DISH for PENETROMETER	A090P002H
SLOTTED WEIGHT for PENETROMETER - 50 g	A090P003H

SLOTTED WEIGHT for PENETROMETER - 100 g

A090P004H

LUBRICATING GREASE PENETRATION CONE

A090P005H

BITUMEN PENETROMETER - Semi-Automatic

A090S00XU