

KINEMATIC VISCOMETER

Code : A435



- Used to determine the kinematic viscosity of the bitumen and road oil.
- Equipped with electronic thermostat, thermometer and a heating unit.
- Equipped with water circulating unit to provide homogeneous/constant temperature all over the bath.

STANDARDS

ASTM D2170

SUPPLIED WITH

- Viscosity Tube(s) in accordance with ASTM D2170 (To be specified at the time of inquiry).
- Digital Thermostat & Indicator with accuracy of 0.02°C.

SUPPLIED WITH

- Cannon Fenske Opaque Tube for Kinematic Viscometer - ASTM
- Cannon Fenske Tube for Kinematic Viscometer
- Ubbelohde Tube for Kinematic Viscometer

ORDERING INFORMATION

Item Name	Item Code
CANNON FENSKE OPAQUE TUBE for KINEMATIC VISCOMETER - ASTM	A435P001H
CANNON FENSKE TUBE for KINEMATIC VISCOMETER	A435P002H
UBBELOHDE TUBE for KINEMATIC VISCOMETER	A435P003H
KINEMATIC VISCOMETER - ASTM	A435X00AU

OTHER PHOTOS

