

ASPHALT > BITUMEN PROPERTIES

# FLASH and FIRE POINT APPARATUS (Cleveland Open Cup)

Code : A470



- Used to determine the maximum temperature bitumen can be safely heated without the danger of an instantaneous flash in the presence of an open flame.
- Supplied with an electronic heater, a thermoregulator, two thermometers, and a sample cup.

## STANDARDS

ASTM D117 • ASTM D92 • AASHTO T-48

## SUPPLIED WITH

- Electric heater

- Thermoregulator
- Two Thermometers
- Sample Cup

**ORDERING INFORMATION**

| Item Name                                                     | Item Code |
|---------------------------------------------------------------|-----------|
| AUTOMATIC FLASH and FIRE POINT APPARATUS (Cleveland Open Cup) | A470A00XU |
| FLASH and FIRE POINT APPARATUS (Cleveland Open Cup)           | A470H00XU |