


RING AND BALL APPARATUS (Automatic)

Code : A472A



will be added soon ...

TECHNICAL SPECIFICATIONS

will be added soon ...

ORDERING INFORMATION

Item Name	Item Code
AUTOMATIC RING AND BALL APPARATUS	A472A00XU
RING and BALL SOFTENING POINT APPARATUS - ASTM	A472H00AU
RING and BALL SOFTENING POINT APPARATUS - EN	A472H00EU
STEEL BALL for SOFTENING POINT APPARATUS	A472P001H
BALL CENTERING GUIDE for SOFTENING POINT APPARATUS	A472P002H

BRASS RING for SOFTENING POINT APPARATUS

A472P003H