


COMPRESSOMETER/EXTENSOMETER for CONCRETE SPECIMENS

Code : B002P



will be added soon ...

TECHNICAL SPECIFICATIONS

will be added soon ...

ORDERING INFORMATION

Item Name	Item Code
FLEXURAL TESTING MACHINE (C-Type) - 200 kN	B002A20XU
FLEXURAL TESTING MACHINE (C-Type) - 300 kN	B002A30XU
POWER PACK (for Compression Frame) - 400 Bar	B002P000U

POWER PACK (for Flexure Frame) - 400 Bar	B002P001U
POWER PACK (for Compression Frame) - 400 Bar (Without Controlling Unit)	B002P002U
POWER PACK (Servo-Motor, for Compression Frame) - 400 Bar	B002P003U
HYDRAULIC UNIT (for Compression Frame) - 400 Bar	B002P004U
POWER PACK (for Semi-Automatic Compression Frame) - 400 Bar	B002P005U
POWER PACK with PILOT VALVE - 400 Bar (Without Controlling Unit)	B002P006U
UPGRADE PACK for MODULUS OF ELASTICITY MEASUREMENT	B002P007H
POWER PACK with PILOT VALVE - 400 Bar	B002P008U
FLEXURAL FRAME (C-Type) - 200 kN	B002P20XH
FLEXURAL FRAME (C-Type) - 300 kN	B002P30XH