

## CAPPING SET

Code : B015



- Comprises a melting pot and a capping apparatus.
- The melting pot is designed to melt capping compounds and features an adjustable temperature control via a thermoregulator.
- A built-in pilot lamp indicates the heating status.
- The pot is constructed with a double-walled frame for enhanced durability and insulation.

### STANDARDS

ASTM C31 • ASTM C192 • EN 12390-3 • AASHTO T23

### SUPPLIED WITH

- Capping apparatus (Used for Ø 100 mm and Ø 150 mm cylinders)

## **TECHNICAL SPECIFICATIONS**

- Melting pot (4.5 lt):
  - Thermostatically Controlled
  - Insulated double-wall
- 220 - 240 V / 50 - 60 Hz

## **ORDERING INFORMATION**

| Item Name                                 | Item Code |
|-------------------------------------------|-----------|
| CAPPING FRAME (for Ø 100 mm and Ø 150 mm) | B015P001H |
| MELTING POT - 5 lt                        | B015P002U |
| MELTING POT - 10 lt                       | B015P003U |
| CAPPING FRAME (for Ø 110 mm and Ø 160 mm) | B015P004H |
| CYLINDER CAPPING SET                      | B015X00XU |

## **OTHER PHOTOS**

