

REBAR DETECTOR

Code : B050



- Used to detect the most frequent materials found behind walls, ceilings and floors, including ferrous and non-ferrous metals, AC live wires and studs.

TECHNICAL SPECIFICATIONS

Specification	Value
Detectable materials	Electricity, ferrous and non-ferrous metals, wood and metal substructures in drywalls
Detection depth, steel, max.	100 mm
Detection depth, copper, max.	80 mm
Detection depth of live cables, max.	50 mm
Detection depth of wood, max.	25 mm

ORDERING INFORMATION

Item Name	Item Code
REBAR DETECTOR (Bosch)	B050X001B
REBAR DETECTOR (Proceq, Profoscope+)	B050X002B
COVER METER (Proceq, Profometer 650 AI)	B050X003B