

POKER VIBRATOR

Code : B128



- The lab-type poker vibrator is used for internal compaction for the concrete specimens.
- The flexible shaft is 1 meter long.

TECHNICAL SPECIFICATIONS

- Ø 25 mm Flexible shaft (1 m long)
- Direct Connection to the Motor
- 220 - 240 V / 50 - 60 Hz

ORDERING INFORMATION

Item Name	Item Code
POKER VIBRATOR (Battery-Operated) - Ø 25 mm	B128X25XB

POKER VIBRATOR - Ø 25 mm	B128X25XU
POKER VIBRATOR - Ø 28 mm	B128X28XU
POKER VIBRATOR - Ø 35 mm	B128X35XU
POKER VIBRATOR - Ø 45 mm	B128X45XU