

CEMENT > FINENESS OF CEMENT

BLAINE APPARATUS (Automatic)

Code : C035A



will be added soon ...

TECHNICAL SPECIFICATIONS

will be added soon ...

ORDERING INFORMATION

Item Name	Item Code
AUTOMATIC BLAINE AIR-PERMEABILITY APPARATUS - Digital	C035A00XU
BLAINE AIR-PERMEABILITY APPARATUS - Manual	C035H00XH
FILTER PAPER for C035	C035P001H
GLASS U-TUBE MANOMETER for C035	C035P002H
SEMI-AUTOMATIC BLAINE AIR-PERMEABILITY APPARATUS - Digital	C035S00XU