

CEMENT > WORKABILITY & CONSISTENCY

VICAT APPARATUS (Automatic)

Code : C090A



will be added soon ...

TECHNICAL SPECIFICATIONS

will be added soon ...

ORDERING INFORMATION

Item Name	Item Code
AUTOMATIC VICAT APPARATUS - ASTM	C090A00AH
AUTOMATIC VICAT APPARATUS - EN	C090A00EU
VICAT APPARATUS - ASTM	C090H00AH
VICAT APPARATUS - EN	C090H00EH
VICAT MOULD - ASTM	C090P001H

VICAT MOULD - EN	C090P002H
VICAT INITIAL SETTING TIME NEEDLE - ASTM	C090P003H
VICAT INITIAL SETTING TIME NEEDLE - EN	C090P004H
VICAT FINAL SETTING TIME NEEDLE - EN	C090P006H
CONICAL PENETRATION NEEDLE - EN 13279-2	C090P007H
VICAT WEIGHT - 700 g	C090P008H
CONSISTENCY PLUNGER for VICAT APPARATUS - Ø 10 mm	C090P009H
GRADUATED MEASURING PLATE for VICAT APPARATUS	C090P010H