

COMPRESSOR

Code : G701



- Compressors in engineering testing laboratories are essential for providing a reliable source of high-pressure air. They power pneumatic tools, automate testing equipment, assist in sample preparation (e.g., cleaning or drying), and ensure consistent operation of systems requiring compressed air. This improves efficiency, accuracy, and repeatability in various testing procedures.

ORDERING INFORMATION

Item Name	Item Code
COMPRESSOR - 8 bar / 25 lt	G701X025U
COMPRESSOR - 8 bar / 50 lt	G701X050U
COMPRESSOR - 10 bar / 300 lt	G701X300V