

PH METER

Code : G705



- A pH meter in engineering testing laboratories is used to measure the acidity or alkalinity of materials, especially in soil and water testing. It helps assess chemical properties that affect material behavior, ensuring compliance with standards and supporting accurate analysis for construction and environmental decisions.

ORDERING INFORMATION

Item Name	Item Code
PH METER (Bench Type)	G705X001U
PH INDICATION PAPER	G705X002H
PH METER (Pocket Type)	G705X003B
PH METER (with Probe)	G705X004B

PH CALIBRATION SOLUTION (pH 4 / 7 / 10)

G705X004H