

MIXING BOWL

Code : G722



- A mixing bowl in engineering testing laboratories is used for preparing and mixing materials like soil, cement, or mortar. It provides a clean, controlled space for uniform mixing, ensuring consistent and accurate test samples.

ORDERING INFORMATION

Item Name	Item Code
MIXING BOWL with HEATER - 5 lt	G722X005U
MIXING BOWL - Ø 12 cm	G722X012H
MIXING BOWL - Ø 20 cm	G722X020H
MIXING BOWL - Ø 28 cm	G722X028H