

## SPIRIT LEVEL

Code : G764



- A spirit level in engineering testing laboratories is used to check the horizontal or vertical alignment of surfaces or equipment. It ensures accurate positioning during sample preparation, equipment setup, and testing procedures, improving the reliability of results.

### ORDERING INFORMATION

| Item Name             | Item Code |
|-----------------------|-----------|
| SPIRIT LEVEL - 20 cm  | G764X020H |
| SPIRIT LEVEL - 30 cm  | G764X030H |
| SPIRIT LEVEL - 60 cm  | G764X060H |
| SPIRIT LEVEL - 80 cm  | G764X080H |
| SPIRIT LEVEL - 120 cm | G764X120H |

SPIRIT LEVEL - 180 cm

G764X180H