

## QUANTAB CHLORIDE TITRATOR STRIPS (Type 1175)

Code : G815



- The Quantab Chloride Titrator 1175 type in engineering testing laboratories is used for determining the chloride content in water, soil, or concrete samples. It is a quick and efficient method for assessing the chloride levels, which is crucial for evaluating corrosion potential, especially in materials like reinforced concrete.

### **ORDERING INFORMATION**

| Item Name                           | Item Code |
|-------------------------------------|-----------|
| QUANTAB CHLORIDE TITRATOR 1175 TYPE | G815X00XH |