

SHATTER INDICES TESTING MACHINE

Code : M015



- Used to determine the resistance to breakage after dropping the sample with specified conditions mentioned in standards
- The box has a lifting mechanism by pulleys and 2 door hinged lengthwise with a latch for rapid opening
- Rigid base plate with 4 sides preventing loss of coke

TECHNICAL SPECIFICATIONS

- Equipped with pulleys and hinged doors
- Rigid base plate with 4 sides preventing loss of coke

ORDERING INFORMATION

Item Name	Item Code
SHATTER INDICES TESTER	M015X00XH