

GRADING and WEIGHING SYSTEM

Code : M041



- Used wherever accepting the ore to the factory
- There are 4 functions in the system:
 - **Lifting the sample:**
 - Capacity of the crane winch is 2000 kg
 - Travel beam of the winch is limited by switchers
 - Manual controlling device of the winch to be used by operator
 - **Feeding:**
 - The bigbag is filled to the feeding tank by using the crane
 - The feeding tank is carried on to the sieving system
 - The bottom door is opened hydraulically by the operator
 - **Grading:**
 - Grading is done by the 3 deck screening

- Each size is collected in the collector

- **Weighing:**

- Different sizes of the materials are weighed by the platform scale

TECHNICAL SPECIFICATIONS

- **The system is complete with:**

- Platform scale with 2000 kg capacity (digital)
- Crane winch and supporting system with control unit
- Three deck grading machine with each deck 1000 x 2000 mm size
- Feeding tank with capacity of 1500 kg with hydraulic bottom door
- Trolleys made of steel construction (4 pcs)

ORDERING INFORMATION

Item Name	Item Code
GRADING and WEIGHING SYSTEM	M041X00XV