

PUNCH MACHINE

Code : M043



- Used to take samples from metal sheets up to 2 mm thick
- Hydraulic system with capacity of 1000 kN
- Supplied with two moulds that can take samples with diameters of 57.3 mm and 64.5 mm
- The system has a safety valve for overload
- Equipped with safety guard door that stops the machine once opened

TECHNICAL SPECIFICATIONS

- Capacity: 1000 kN
- Equipped with safety guard door
- 380 V (3 phase)

ORDERING INFORMATION

Item Name	Item Code
PUNCH MACHINE	M043X001U
STEEL SHEET SAMPLING MACHINE	M043X00XU