

ALUMINUM SAMPLING MACHINE

Code : M062



- Used to prepare circular samples from aluminum ingots for XRF analysis
- Fully automatic, with 3 functions:
 - Drilling
 - Polishing
 - Cutting
- Ingot is placed on the support system and all the a/m functions are done automatically in a short time by an operator
- The sample taken has a desired diameter and thickness with smooth surface

TECHNICAL SPECIFICATIONS

- Fully automatic, with 3 functions:
 - Drilling

-
- Polishing
 - Cutting
 - Ingot is placed on the support system and all the a/m functions are done automatically in a short time by an operator

ORDERING INFORMATION

Item Name	Item Code
ALUMINUM SAMPLING MACHINE	M062X00XV