

# CYLINDER CRUSHER

Code : M415



- Made of hardened steel cylinders. The distance between the cylinders is adjustable.
- The materials feed size is < 20 mm, and it is supplied with overload protection with 2 springs.
- Supplied with feed hopper and collector.
- Supplied complete with the electrical motor.

## **TECHNICAL SPECIFICATIONS**

- Material feed size < 20 mm
- Cylinders are steel and hardened.
- Distance between cylinders is adjustable
- Overload protection with 2 springs
- Feed hopper and collector
- 380 V / 50 Hz

**ORDERING INFORMATION**

| Item Name                | Item Code |
|--------------------------|-----------|
| SPARE CYLINDER for G-415 | M415P001H |
| CYLINDRICAL CRUSHER      | M415X00XU |