

SOIL > MOISTURE-DENSITY RELATIONSHIP

SOIL COMPACTOR - EN (Automatic)

Code : T055 (EN)



will be added soon ...

TECHNICAL SPECIFICATIONS

will be added soon ...

ORDERING INFORMATION

Item Name	Item Code
AUTOMATIC SOIL COMPACTOR - ASTM [60 Hz]	T055A00AK
AUTOMATIC SOIL COMPACTOR - ASTM [50 Hz]	T055A00AT
AUTOMATIC SOIL COMPACTOR - EN	T055A00EU
NOISE-REDUCTION and SAFETY CABINET for T-055/ASTM	T055P001H
TURNABLE BASE with STEPPER MOTOR for T055 (ASTM Model)	T055P002H

RUBBER WEDGE - Complete for T055 (EN Model)	T055P003H
SPARE GEAR BOX for T055 (ASTM Model)	T055P004H
SPARE HAMMER for T055 (ASTM Model)	T055P005H
ELECTRONIC CONTROLLING UNIT for T055 (EN Model)	T055P006H
SPARE HAMMER SPRING for T055 (EN Model)	T055P007H
TURNABLE BASE with STEPPER MOTOR for T055 (EN Model)	T055P008H
ELECTRONIC BOARD (PCB) for T055 - ASTM	T055P009H
ASTM STANDARD RAMMER for T055A00EU	T055P010H