

SOIL > MOISTURE DETERMINATION

UNIVERSAL CARBIDE TESTER (Speedy Moisture)

Code : T068



- Used to determine the moisture in the soil, sand and aggregate particles. The machine gives quick and accurate results. The tester is suitable for both lab and field work.
- The mechanism carried out by the machine is through measuring the gas emission from the reaction between water and calcium carbide.
- The tester is equipped with a portable electronic balance, steel balls and a carrying case.

SET

- Calibrated Pressure Bottle
- Moisture Meter
- Precise Digital Scale
- Spherical balls
- Sample Cups (2 ea)

- Hammer
- Chisel
- Sampling Spoon
- Sample Grinding Bowl
- Set of 25 ampoules of Calcium Carbide
- Digital Timer/Stopwatch
- Calibration set
- Cleaning Brush
- Carrying Case

ORDERING INFORMATION

Item Name	Item Code
TEST POWDER (Pack of 25)	T068P001H
TEST POWDER (Pack of 100)	T068P002H
UNIVERSAL CARBIDE TESTER (Speedy Moisture Tester)	T068X00XB