

WATER LEVEL INDICATOR

Code : T077



To be added soon ...

TECHNICAL SPECIFICATIONS

To be added soon ...

ORDERING INFORMATION

Item Name	Item Code
WATER LEVEL INDICATOR - 50 m	T077X050B
WATER LEVEL INDICATOR - 100 m	T077X100B
WATER LEVEL INDICATOR - 150 m	T077X150B
WATER LEVEL INDICATOR - 200 m	T077X200B
WATER LEVEL INDICATOR - 300 m	T077X300B