


WEIGHT SET for CONSOLIDATION and DIRECT SHEAR

Code : T215



- Weight set to be used with consolidation (oedometer) testers and direct shear testers.

TECHNICAL SPECIFICATIONS

- Made of galvanized steel

ORDERING INFORMATION

Item Name	Item Code
SLOTTED WEIGHT for OEDOMETER / DIRECT SHEAR - 1 kg	T215P001H
SLOTTED WEIGHT for OEDOMETER / DIRECT SHEAR - 5 kg	T215P005H

SLOTTED WEIGHT for OEDOMETER / DIRECT SHEAR - 10 kg	T215P010H
---	-----------

WEIGHT SET - 50 kg	T215X050H
--------------------	-----------

WEIGHT SET - 64 kg	T215X064H
--------------------	-----------

WEIGHT SET - 80 kg	T215X080H
--------------------	-----------